

Performance Test Report
- AngelGuard GHB Drug Test Sticker
(Detection Limit and Beverage Reactivity)

Approval

☐ Controlled Copy : (No.) ☐ Uncontrolled Copy

Owner:	
Written by:	Sign
Reviewed by:	Sign
Approved by:	Sign

Copyright © 2025 Philmedi All rights reserved.

1. Purpose

This report aims to verify the performance of AngelGuard GHB Drug Test sticker product and to evaluate the validity of its detection capability.

- Key Test Item 1: Detection Limit
- Key Test Item 2: Beverage Reactivity

2. Experimental Design

- The AngelGuard™ GHB drug test sticker (hereinafter “the product”) was evaluated for detection performance using GHB (Gamma-Hydroxybutyric acid) as the standard substance.
- The buffer for diluting GHB consisted of 13 types of beverages. Positive samples were prepared by adding GHB standard solution at 1% concentration; negative samples contained no GHB.
- The detection limit of the product was determined.
- To confirm improvements in detection ink performance, beverage-specific tests were conducted.

During beverage-specific tests, it was also observed whether color changes of positive samples occurred within the waterproof layer.

3. Materials

3.1. GHB Standard Substance

3.2. Buffer Standard Substance (20% ethanol in DW)

3.3. Nail sticker product

3.4. 13 types of beverages (tap water, 20% ethanol in DW, vodka, makgeolli, beer, whiskey, red wine, white wine, orange juice, cola, cider, ion beverage, Red Bull)

Test Report

Performance Test Report
(Detection Limit and Reactivity)

4. Methods

4.1. Storage Conditions

- Temperature: 2 ~ 30 °C (room temperature)
- Environment: Stored inside OPP pouch, outer paper packaging
- Storage place: Philmedi Desiccator (PG-E-2022-07)

4.2. Preparation of Samples

- Gamma-hydroxy-butyric acid sodium salt Standard (GHB CoA attached)



**Certificate of Analysis
Reference Material**

QC-CA-538
538.1B1.2
Page: 1/3

Lipomed Document QC-CA-538
Version: 006-26.Jan.2022


Supersedes: 005-12.Dec.2014

Chemical Name: **GHB-Na**
γ-Hydroxy-butyric acid sodium salt

Lot No: 538.1B1.2
Art. No: GHB-538-NA

Chemical Formula: C₄H₇O₃Na
Molwt: 126.09

CAS Registry No: 502-85-2



TEST	SPECIFICATION	RESULT
1. Appearance	white to off-white crystalline powder	conforms
2. Purity	HPLC > 98.5 % The purity is calculated from the distribution of 8 HPLC analyses with a 95% level of confidence.	99.950 ± 0.050 %
3. Optical Rotation	N/A (no chiral center)	N/A
4. Water Content	Karl Fischer ≤ 1.0 %	0.2 %
5. Identity	IR UV: in methanol No maxima	IR corresponds UV: in methanol No maxima

Storage Conditions: For maximum stability store air-tight at 2 - 8°C in a dark location.

Lipomed certifies that this standard meets the specifications stated in this certificate and warrants this product to meet the stated acceptance criteria through the retest date when stored unopened as recommended.

FOR ANALYTICAL PURPOSES ONLY: NOT FOR HUMAN OR ANIMAL USE!

Issued by Dr. L. Prévot



Date sign: Arlesheim,
January 26, 2022

LIPOMED AG
Fohrmattenweg 4
4144 Arlesheim
Switzerland
Tel. +41 (0) 702 02 00
Fax +41 (0) 702 02 20

LIPOMED GmbH
Hegenheimer Str. 2
79576 Weil am Rhein
Germany
+49 7621 1693 473
+49 7621 1693 474

LIPOMED Inc.
150 Cambridgeport Drive, Suite 705
Cambridge, MA 02140
U.S.A.
+1 (617) 577 7222
+1 (617) 577 1776
www.lipomed.com
lipomed@lipomed.com





GHB Certificate of analysis(CoA)

4.2.1. Preparation of Positive Samples with 13 Beverages

- Dispense 1,000 µL of each beverage and add GHB standard to achieve 1% concentration.

4.3. Method 1 – Detection Limit Test for GHB

- Prepare positive samples with 1% GHB in 20% ethanol.
- Leave for 1 minute, wipe off, and observe changes in waterproof ink region.
- Apply identical positive samples to stickers of different shapes to check consistency of color change.
- Record results by photographs

4.4. Method 2 – Detection of GHB in 13 Beverages

- Detect positive and negative samples according to product instructions.
- Record results by photographs.

4.5. Product Usage Procedure

- 1) Each kit contains 10 test spots (stickers)
- 2) Detach one sticker from the release PET film and attach it to a fingernail.
- 3) If a beverage is suspected to contain GHB, apply it to the sticker using a hand or a straw.
- 4) After 30 seconds (or 1 minute), shake or wipe off the applied sample.
- 5) Determine positive/negative result according to interpretation criteria.

5.	Test Conditions
----	-----------------

5.1. Location: Philmedi Corporate Research Institute (33 Sagimakgol-ro 62beon-gil, Jungwon-gu, Seongnam-si, Gyeonggi-do, Korea)

5.2. Date: September 10, 2025 (1 day)

5.3. Environment

-Temperature 26.3 ± 5 °C, Humidity 46 ± 2 %

5.4. Performed by: Researcher Sejin Kim

6. Interpretation of results

Visual inspection criteria - Naked eye			
Detection schematic		Result	Mark
		Positive	(++) (+)
		Negative	(-)

Test Report



PMD-23-DTSN-DFR-05/R0/2025.09.10

Performance Test Report
(Detection Limit and Reactivity)

7. Results

7.1. Detection Limit Test

No.	Sample Type	Result	Mark
1	20% ethanol 1% GHB (Positive)		++
2	20% ethanol 0% GHB (Negative)		-
- All sticker designs showed identical positive color change.			

Test Report



PMD-23-DTSN-DFR-05/R0/2025.09.10

Performance Test Report
(Detection Limit and Reactivity)

7.2. GHB Detection Test with 13 Beverages

No.	Beverage	Positive Sample	Negative Sample	Positive Reactivity	Negative Reactivity	Waterproof Ink Damage
1	Tap Water			++	-	N/A
2	20% ethanol			++	-	N/A
3	Vodka			++	-	N/A
4	Makgeolli			++	-	N/A
5	Beer			++	-	N/A
6	Whiskey			++	-	N/A

Test Report



PMD-23-DTSN-DFR-05/R0/2025.09.10

Performance Test Report
(Detection Limit and Reactivity)

7	Red Wine			+	-	N/A
8	White Wine			+	-	N/A
9	Orange Juice			-	-	N/A
10	Cola			+	-	N/A
11	Cider			++	-	N/A
12	Ion Beverage			+	-	N/A
13	Red Bull			-	-	N/A

8.	Interpretation
----	----------------

8.1. Detection Limit

- Positive response confirmed at 1% GHB concentration in 20% ethanol
- Stickers of different shapes showed identical color change at detection limit concentration
- Therefore, beverage-specific tests were standardized to one sticker shape.

8.1. Beverage-Specific Detection

- False negatives occurred with orange juice and Red Bull
(False negative = when positive sample was applied but negative result observed)
- Red wine showed coloration in negative samples, but positive/negative distinction was still possible by naked eye.

9.	Conclusion
----	------------

- The product showed positive color change from 1% (10 mg/mL) GHB concentration
- In the 13-beverage GHB detection test, false negatives were observed in orange juice and Red Bull
- All other beverages showed normal positive detection
- Although the intensity of color change varied by beverage, positive/negative results were distinguishable by naked eye